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
GB

Instructions for use

1. Expose the cervix by a speculum and cleanse it with a disinfectant liquid.
2. Perform a hysteroscopy to determine uterus depth.
3. Dilate cervix (5 ½ Hegar)
4. Take the manipulator out of the sterile pouch and check balloon integrity (inflate and deflate with a syringe).
5. Lubricate the tip of the catheter and the balloon; insert the manipulator gently into the uterine cavity through the cervical canal. **Caution:** the catheter must be placed 1 cm underneath detected maximum depth. Push the catheter forward until the adjustable positioner (already placed 1 cm under detected depth) stops against the external cervix.
6. Gently inflate the balloon with 10 cc air to reach a maximum expansion of 20 mm.
7. Gently push the cervical stopper forward against the cervix so that the catheter is secured between the inflated balloon and the cervical stopper. Lock the manipulator with the thumb. The catheter is now placed in the correct position and allows a safe and atraumatic uterus manipulation.
8. Inject the contrast medium by the second way of the catheter to display the uterus and the Fallopian tubes.
9. At the end of the procedure deflate the balloon by connecting the syringe to the one-way valve. Unlock the manipulator with the thumb and remove the catheter gently.

Caution: to be used only after performing a hysteroscopy in order to avoid any uterus perforation. During handling, avoid squeezing the tabs underneath the locking system. Read the instructions before using the device. Do not reuse. Use only if the package is not damaged. Use immediately after opening. Dispose after use according to the laws in force. For use only under trained medical staff.

IF RE-STERILIZED, THE DEVICE COULD GET DEFORMED

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