## PM-CARE®HYSTEROMETRE

HYSTEROMETER



The **PM-CARE® HYSTEROMETRE** hysterometer is used to assess uterine depth during gynaecological examinations, in particular for the insertion or removal of a contraceptive device.

## AVAILABLE REFERENCE(S)

Reference	Diameter (Fr)	Diameter (mm)	Colour	
PS1251000	10	3.3	Blue	
PS1251200	12	4	Grey	
PS1251400	14	4.6	Green	

## **ABOUT THIS MEDICAL DEVICE**

1 <sup>st</sup> release date	1999	
Single use	Yes	
Class	ls	
Sterile	Yes	
Sterilisation method	Ethylene oxide sterilisation	
Contains latex	No	
Lifetime	5 years	
Certification	CE	
Packaging	Box of 25 or 10 units in individual sterile wrapping	
Storage	In its original box in a cool and dry place	

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76

prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.



<sup>®</sup> Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.