




# SURGERY





A1511


## Diamond-Gerjy Proctoscope I A1511

It has never been easier – Wing Plast's Diamond-Gerjy Proctoscope with universal light adapter is made from a clear, environmentally friendly plastic and is suitable for all proctological examinations. The Diamond-Gerjy Proctoscope has a sliding shutter on one side and a window on the opposite side to give a greater view and examination area.

Instruments can be used at all angles without obstructing the view of the operator. The sliding shutter can be opened to the desired length depending on the nature of the treatment. It gives an excellent view and at the same time protects sensitive tissues. Local anesthetic can be injected to increase therapeutic options. The Diamond-Gerjy Proctoscope protects sensitive areas in the anal canal and exposes the desired areas clearly in order to facilitate the treatment of rubber band ligation, injection, suture, pexi and tissue sample etc.

The Diamond-Gerjy Proctoscope is excellent for all proctological surgery.

- Adjustable slide
- Larger examination area
- Better view of the entire anal canal
- Protects sensitive skin and mucosa with the slide during treatments
- Greater flexibility with many instruments
- Ergonomically stable construction



## Proctoscope I 1524/1522/1531/1521/1541

Wing Plast's Proctoscopes with a lamp holder is a disposable instrument manufactured from a clear environmentally friendly plastic. The Proctoscopes are available in sizes: Small, Medium and Large and are ergonomically and functionally designed for the comfort of both user and patient. The lamp holder fits most common adapters.

- Stable construction
- Easy to manoeuvre
- Lamp holder for fiber optic adapter
- Optimal illumination in the Proctoscopes distal tip
- Smooth rounded edges
- A smooth transition between the tip of the obturator and the end of the Proctoscope to maximize the comfort for the patient



**H1002**

## Multi-Ligator I H1002

The Multi-Ligator PLUS is a disposable device for the band ligation of haemorrhoids. It is preloaded with 4 latex-free bands, has a transparent nozzle and a trigger guideline. It has never been easier - Place the tip of the Multi-Ligator on the haemorrhoid and gently squeeze the trigger to the marked line, which induces the required amount of suction. Continue to pull the trigger to apply the band. When the trigger is released, the Multi-Ligator is automatically reloaded, ready for immediate use.

- Disposable
- Easy to use, can be used four times before being reloaded
- Treatment only takes 5 minutes
- Can be used 4 times without reloading
- Reduced pain to the patient
- New transparent nozzle
- Better visibility
- Stronger elastic bands
- Cost-effective




**A430**

## White Balance Cup I A430

Wing Plastic's sterile disposable white balance cup is used for white balance of surgeon cameras. It is made of environmentally friendly plastic.

- Stable and user friendly
- Environmentally friendly plastic
- Sterile
- Fast white balancing with high quality
- Saves time for the user
- Protects the optical lens from ambient light in a good way.



1542

## Anoskop I 1542

Wing Plast's Anoskop 1542 for single use is made of environmentally friendly plastic. The anoscope is used in colorectal examinations and is constructed of a material that conducts light very well and provides an optimal view for the examiner. The anoscope's carefully rounded edges and surfaces contribute to a gentle examination for the patient.

- The beveled end provides better access to hemorrhoid treatments.
- Fits with the same light fixture as for rectoscopes and tube proctoscopes.
- Fits with several of the lighting systems on the market today.



1551

## Tube Proctoscope I 1551

Wing Plast's disposable Tube Proctoscope, is manufactured in a clear environmentally friendly plastic with excellent light transmission which gives an optimal, distal, annular light view for the practitioner. The Tube Proctoscope has smooth edges that contribute to a gentle examination for the patient.

- Contoured design for ease of use
- Clearly marked graduated scales
- Fits most popular light heads
- A smooth transition between the tip of the obturator and the end of the Proctoscope to ensure maximum patient comfort
- Exceptional distal tip illumination gives a bright and clear view



1561

1561-1

## Sigmoidoscope I 1561/1561-1

Wing Plastics Rectoscope 1561-1 for disposable uses are manufactured from a environmentally friendly plastic, the tube is constructed from a material which gives an optimal, distal, annular light view for the practitioner. Rectoscope's rounded edges and surfaces contribute to a gentle examination for the patient.

- Contoured finger grip for ease of use
- Clearly marked graduated scales
- Fits most popular light heads
- A smooth transition between the tip of the obturator and the end of the Sigmoidoscope ensures maximum patient comfort
- Exceptional distal tip illumination gives a bright and clear view



1561-PL

## Rectoscope kit with light disposable I 1561-PL

Wing Plastics Rectoscope 1561-PL for Disposable use manufactured from environmental friendly plastic, the tube is constructed of a material which conducts light very well and provides an optimal view for the examiner. Rectoscope's rounded edges and surfaces contribute to a gentle examinationfor the patient.

- Stable and grip-friendly construction
- Clear scale
- Fits different light heads
- Very small transition between obturator and tube
- The light transmission gives the correct tissue color and free visibility





1562

## Insufflation Bellow I 1562

Wing Plast's reusable Insufflation Bellow is an easy-to-use bellow with a expandable ballon. Used for rectoscopy and is provided with a safety net.

Also available as single bellow, art. no. 1566

- High quality
- Suitable for different Light heads
- Easy to hold



1563

## Double Balloon Rekto Disposable I 1563

Wing Plast's Double Bellow, disposable with a safety net for Rectoscopy examinations has the following benefits:

- Time effective
- No cleaning needed
- Better hygiene, reduced risk of spreading infections
- Fits most Light heads
- Latex free



1565

## Rectoscope filter I 1565

Wing Plast's disposable Filter 1565, is a bacteria and vent filter for Sigmoidoscopes. It is placed between the Sigmoidoscope Light Head and the Insufflation Bellow, in order to prevent liquid and other particles from entering into the Insufflation Bellow.

- Prevention of the bacteria and spread of infections to 99.999 % of cases
- Better hygiene for both the examiner and the patient
- Includes 2 Luer Connectors in a 25-pack



1565

## Luer coupling I 1569

Wing Plasts Luer Coupling 1569 is used as a link between a rectoscope filter and a double balloon, make it easier applying the filter.

- Luer lid bracket for smooth attachment of filter 1565
- Reusable
- Easy grip





**A410**

## Forceps I A410

Wing Plast's Tweezers A410 is a practical one and very useful instrument in many different treatment situations. It's easy and efficient to work with, for example when removing sutures, resurfacing and catheter management.

- Stable construction
- Environmentally friendly
- Easy grip
- Sterile



**1550**

## Anal Dilator I 1550

Wing Plast's Anal Dilator for self-treatment, is a dilator for prevention of strictures following proctological surgery and other painful anal conditions. Use of the Anal Dilator makes the anal sphincter muscle relax which reduces pain and prevents spasms in the anal area.

- Made of environmentally friendly plastic
- Prevents tearing and fissures in the anal canal
- No abnormal widening of the sphincter
- Reusable (for the same person)



1527

1528

1529

1528H

## Fiber Light Adapter 1527/1528/1529/1528H

Wing Plast Fiber Light Adapter fits our Proctoscope and gives a good reflex-free light.

Available for other lighting systems

- Good light conductivity
- Standardized fitting
- Easy to manoeuvre
- Stainless steel
- Autoclavable



1528E

1529E

## Fiberljusadapter I 1528E/1529E

Wing Plast Light adapter disposable for Proktoscope gives a very bright light while eliminating cleaning. Advantages against multiple light adapters always provide optimal light and remove the risk of disappearing.

- Prevents bacteria and spread of infection
- Better hygiene for both examiner and patient
- Save time
- Better light (than multi-fiber fiber light adapter)
- Disposable



1564



1564H



1564WA



1564X-VRID

## Rectoscope Light Head

### 1564/1564H/1564WA/ 1564X-TURN

Wing Plastic Rectoscope Light Head gives a sharp reflex-free light. The light head has an openable door and a pin for attaching Double Balloon.

1564X-TURN has a magnifying glass in the lid.

Available with socket for other lighting systems

1564: Storz

1564X-TURN: Storz

1564H: Heine

1564WA: Welch Allyn

- Stainless
- Autoclavable
- Easy to clean



1570

## Lightsource LED I 1570

Wing Plast's LED light source is an incredibly strong light source for its size. It produces a stronger light than a 150W halogen light source and has one life span of about 50,000 working hours.

Now also available as a complete kit mounted on a stand.

- Compact, efficient and economical light source
- Light: 50 Watt LED
- Lifetime of the LED lamp: 50,000 working hours
- Color temperature: 5600 K
- Fits standard adapters
- 3 years warranty
- CE marked
- Approved according to IEC 60601-1



1575

## LED Light Handle 3,7 V I 1575

Wing plastic's light handle is an incredibly powerful light source for its size. It gives a lot clear light. The light handle comes with two batteries and a charger where you can charge both batteries at the same time.

- Compact, efficient and economical light handle.
- Portable
- Lamp warranty 50,000 hours
- Fits adapters with thread
- Battery life: 2.5 - 12 hours (Level 1 - Level 4)
- 1 year warranty



1576

## Handle Rectoscope I 1576

Wing Plastic handle for rectoscope is an excellent tool if you want a better grip studies. Place on the fiber optic cable (1571/1571-L) and tighten using screw function.



1583



1585



1593



1595

## Projector Adapter/ Instrument Adapter

Wing Plast provides adapters for several different sockets on the market. The most common examples can be found below:

1583 - Projector adapter Olympus

1585 - Storz Projector Adapter

1593 - Olympus Instrument Adapter

1595 - Instrument adapter Storz



1571

## Fiberoptic cable | 1571/1571-L

Wing Plast's Fiber light cable 1571 for halogen/ LED lightsource is lightweight, durable and has a very good light conductivity. The cable is flexible and bendable with an internal protection for maintaining performance.

At both ends, there is a thread for the application of suitable adapter to the respective light source/instrument.

- Very good light conductivity
- Can be ordered in a variety of lengths and colors
- Flexible
- Easy to keep clean
- Autoclaveable 134 °C
- Ø 4,8 mm



4120

## Lubricating Gel I 4120

A high quality Lubricating gel for professional use. Thanks to its ingredients it is best suited for Endoscopy, Gastroscopy, Proctoscopy and Gynecology. Available in 82 g and 260 g.

- Has no negative effect on surgical instruments, rubber or plastics.
- Safe from allergies
- Non-irritating Contains no toxic
- Contains no
- Odorless
- Completely water soluble



1910/1920

## Smart Cover I 1910/1920

Smart Cover is a pair of blue boxers for use when performing an examination. These boxers protects the privacy of the patient whilst enabling the doctor/nurse the necessary working area. Made of non-woven material which is classified as environmental friendly.

- Stretch seams
- Easy to take on and off
- Unisex
- Perforations from front to back
- Suitable for both Gynecological and Colorectal examinations



Since 1987 our Proctoscope have been considered easy to use and stably constructed with a smooth transition between the tip of the obturator and the end of the Proctoscope. This ensures maximum comfort for the patient. Our latest alternative, Dimond-Gerjy Proctoscope is a brilliant solution for healthcare professionals, with an adjustable slide it gives an excellent view and greater examination area. This provides a better view through the anal canal, it protects sensitive skin and mucosa and provides flexibility with many instruments during treatment. It has never been easier. Additionally, Wing Plast has a superb solution for a quicker and easier way to kill haemorrhoids. The Multi-Ligator! Disposable, ready to use, cost effective and reduced pain for the patient, why not give it a try? You can find much more in our catalogue, Tweezers, Sigmoidoscopies, and Filters etc. We may just have the product you are looking for and we will be happy to help you!

**If you have an innovative idea in healthcare and want help with development or manufacturing, please contact us, we'll help you get your idea into reality.**



**Wing Plast** - *Quality Medical Instruments*

Phone +46 (0)371 330 30, Fax +46 (0)371 330 31 , e-mail: [info@wingplast.se](mailto:info@wingplast.se), [www.wingplast.se](http://www.wingplast.se)

RMP MED Steffen Roßberg  
Thomas-Mann-Straße 16

01219 Dresden

# Bestellfax 0351 4188 2297

Telefon  
0351 4188 2291

Fax  
0351 4188 2297

Mail  
info@rmp-med.shop

Online  
www.rmp-med.shop

Artikel-Nr.	Bezeichnung	Größe / Variante	Menge

Ihre Anschrift

ggf. abweichende Lieferanschrift

Allgemeine Geschäftsbedingungen: Lieferung auf Rechnung zu angegebenen Konditionen, zuzüglich anteiliger Versandkosten 6,30 € (7,50 € brutto) pro Bestellung. Zahlungsbedingungen: 14 Tage ohne Abzug. Angebot freibleibend! (vollständige AGBs unter [www.rmp-med.shop](http://www.rmp-med.shop), Stand August 2022)

RMP MED Steffen Roßberg medizinische Produkte e.K.

Rechnungsanschrift:  
Thomas-Mann-Straße 16  
01219 Dresden  
DEUTSCHLAND

Telefon: +49 351 4188 2291  
Fax: +49 351 4188 2297  
E-Mail: [info@rmp-med.shop](mailto:info@rmp-med.shop)  
Web: [www.rmp-med.shop](http://www.rmp-med.shop)

Registergericht Dresden  
HRA 6587  
Ust.-ID-Nr.: DE253156735  
EORI-Nr. DE785413139402371

Ostsächsische Sparkasse Dresden  
DE81 8505 0300 0221 0994 76  
OSDDDE81XXX