

# Bold

Aurelia® Bold® Black Nitrile  
Powder Free Examination Gloves



**Durable Black Nitrile  
Maximum Strength**

- 5 mil thickness for improved puncture resistance
- Fully-textured finish for improved gripping action
- 100% latex free, black color



**AURELIA®**  
Sign of Comfort™

[www.aureliagloves.com](http://www.aureliagloves.com) | 1-877-AURELIA (1-877-287-3542)

# Bold

## Product Description

Aurelia® Bold® Black Nitrile (Non-Latex) Powder Free Examination Gloves are made from 100% Nitrile (Non-Latex), ambidextrous, black color, fully textured and non-sterile. Aurelia gloves conform and are tested to the highest international standards. For more information visit [www.aureliagloves.com](http://www.aureliagloves.com)



## Aurelia® Bold® Specifications

### TENSILE TEST RESULTS (unaged)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	26.1	550
Small	26.1	550
Medium	26.5	550
Large	26.1	595
Extra Large	26.2	580
XX-Large	26.2	580
ASTM D6319-00	14 MPa minimum	500% minimum

### AGED TENSILE TEST RESULTS (7 Days at 70° C)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	24.0	550
Small	24.0	550
Medium	24.0	555
Large	24.0	585
Extra Large	23.7	570
Extra Large	23.7	570
ASTM D6319-00	14 MPa minimum	500% minimum

## DIMENSIONS

	Length (inches/mm)	Thickness (Single-wall)	
		Palm (mils/mm)	Finger (mils/mm)
<b>Extra Small</b> Aurelia Bold ASTM D6319-00	9.4 / 240 8.7 / 220 minimum	5.0 / 0.13 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
<b>Small</b> Aurelia Bold ASTM D6319-00	9.4 / 240 8.7 / 220 minimum	5.0 / 0.13 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
<b>Medium</b> Aurelia Bold ASTM D6319-00	9.4 / 240 9.1 / 230 minimum	5.0 / 0.13 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
<b>Large</b> Aurelia Bold ASTM D6319-00	9.4 / 240 9.1 / 230 minimum	5.0 / 0.13 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
<b>Extra Large</b> Aurelia Bold ASTM D6319-00	9.4 / 240 9.1 / 230 minimum	5.0 / 0.13 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
<b>XX-Large</b> Aurelia Bold ASTM D6319-00	9.4 / 240 9.1 / 230 minimum	5.0 / 0.13 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum

## Ordering Information

To order from your preferred distributor, specify size and item number. Packing: 100 pieces per dispenser box, 10 dispenser boxes per shipping carton.

Extra Small	Small	Medium	Large	Extra Large	XX-Large
73995	73996	73997	73998	73999	73990

Authorized Aurelia Distributor:



AURELIA Div. Supermax Inc.

1899 Sequoia Drive, Aurora, IL 60506 | Toll Free: 1-877-AURELIA | Tel: 1-630-898-8886 | Fax: 1-630-898-8855  
 Web: [www.aureliagloves.com](http://www.aureliagloves.com) | Email: [info@aureliagloves.com](mailto:info@aureliagloves.com) | Facebook: [www.facebook.com/aureliagloves](http://www.facebook.com/aureliagloves)