



AURELIA® Refresh®

Examination Gloves



99225	(XS)	ø 100	10 x 100
99226	(S)	ø 100	10 x 100
99227	(M)	ø 100	10 x 100
99228	(L)	ø 100	10 x 100
99229	(XL)	ø 100	10 x 100



Latex



Powder Free



Green



6.0 +/- 0.2 g



Honeycomb Textured




3 Years

Premium Natural Latex examination gloves that offer maximum strength, comfort, dexterity and sensitivity.

- Polymer Coated Powder-Free Latex
- Peppermint Scented
- Non-Sterile
- Beaded cuff
- Ambidextrous
- Honeycomb Textured
- AQL 1.5

Sizes & Dimensions



XS	70 ± 10 mm	min. 240 mm
S	80 ± 10 mm	min. 240 mm
M	95 ± 10 mm	min. 240 mm
L	110 ± 10 mm	min. 240 mm
XL	≥ 110 mm	min. 240 mm

Single wall thickness



- ① min. 0.15 mm
- ② min. 0.12 mm
- ③ min. 0.08 mm

Physical properties

	Unaged	Aged
Tensile Strength	min. 18.0 MPa	min. 14.0 MPa
Elongation	min. 650 %	min. 500 %
Strength at Break	min. 6N	min. 6N

Quality Signs & Certificates

- 89/686/EEC, 93/42/EWG
- ISO 9001:2008, ISO 13485, CMCAS 13485 (SGS)
- EN 455 1-4
- EN 420 : 2003 + A 1 : 2009
- EN 374 2-3, EN 388
- EN 1186: Food Contact Approved (SGS)

