

AURELIA® Bold MAX®

Examination Gloves



9789A7 (M)	ø 50	10 x 50
9789A8 (L)	ø 50	10 x 50
9789A9 (XL)	ø 50	10 x 50
9789A0 (XXL)	ø 50	10 x 50

Premium quality black nitrile examination gloves with an excellent diamond grip finish offering a greater tensile strength and increased protection than both Latex and Vinyl

- 6 Mil thickness
- 7.5g in weight
- Diamond Textured
- Latex Free
- Ambidextrous
- AQL 1.5



Nitrile



Powder Free



Black



7.5 +/-0.2g




Diamond-Textured



5 Years

Sizes & Dimensions



M	95 ± 10 mm	min. 240 mm
L	110 ± 10 mm	min. 240 mm
XL	120 ± 10 mm	min. 240 mm
XXL	≥ 120 mm	min. 240 mm

Single wall thickness



- ① min. 0.20 mm
- ② min. 0.16 mm
- ③ min. 0.08 mm

Physical properties

	Unaged	Aged
Tensile Strength	min. 14.0 MPa	min. 14.0 MPa
Elongation	min. 500 %	min. 400 %
Strength at Break	min. 6N	min. 6N

Quality Signs & Certificates

- PPE Regulation (EU) 2016/425, 93/42/EWG
- ISO 9001:2008, ISO 13485, CMDCAS 13485 (SGS)
- EN 455 1-4
- EN 420 : 2003 + A 1 : 2009
- EN 374 1-5:2016, EN388
- EN 1186: Food Contact Approved (SGS)

